UNICSKIN

YOUR EFFECTIVE SOLUTION

UNIC**QUANTUM-V GuaSha™ 5.0**

5 TECHNOLOGIES · 3 SMART MODES

DRAIN · FIRM · GLOW



BIPOLAR RADIOFREQUENCY | **EMS TECHNOLOGY** | **RED LED light** | **GREEN LED light** | **NIR LIGHT**

ILLUMINATES & UNIFIES THE TONE | SMOOTHES LINES & WRINKLES | DEFINES THE FACIAL CONTOUR | FIRMING & TIGHTENING | DEPUFFS & REDUCES SWELLING

USER MANUAL

LET'S BEGIN

Congratulations! You have taken the first step towards firmer, smoother, and visibly rejuvenated neck skin, with UNICSKIN, global experts in medical beauty technology and at-home skincare.

Before you start using this device, carefully read all instructions and precautions and use it only for the purpose described in this manual.

1 DEVICE, 5 TECHNOLOGIES, 3 MODES, 5 BENEFITS

UNIC**QUANTUM-V GuaSha™ 5.0** reinvents the ancient Gua Sha ritual with cutting-edge technology. It combines five technologies in a single device: Bipolar Radiofrequency, Low-Frequency EMS, Red LED, Green LED, and Near-Infrared Light, to improve firmness, contour, tone, and radiance of the face and neck.

Compact, lightweight, and 100% wireless, it offers a professional-inspired facial ritual in just 6 minutes a day. Its three intelligent modes work progressively: DRAIN decongests and prepares, FIRM firms and defines contours, and GLOW illuminates and evens skin tone, for visibly firmer, smoother, and more radiant skin.

MAIN BENEFITS

- FIRMING & TIGHTENING
- IDRAINS & REDUCES SWELLING
- BRIGHTENS & EVENS SKIN TONE
- DEFINES THE FACIAL CONTOUR
- REDUCES THE APPEARANCE OF BLEMISHES
- SMOOTHS LINES & WRINKLES

THE SCIENCE BEHIND FIRMER, MORE RADIANT SKIN

This device combines bipolar radiofrequency, low-frequency EMS, and red, green, and NIR LED light to treat the face and neck comprehensively.

Each technology works in a complementary way to help treat the main visible signs of aging in the area: loss of firmness, wrinkles, expression lines, swelling, uneven tone, and spots.

Each technology has a specific function:

Bipolar Radiofrequency (RF)

It emits deep, controlled heat that raises the temperature of the dermis to stimulate fibroblasts and collagen synthesis, helping to improve the appearance of firmness, elasticity, and facial contour. Integrated into FIRM mode.

Low Frequency EMS

Gentle electrical stimulation helps tone and define facial muscles, promoting a more lifted appearance and helping to reduce the appearance of puffiness. It is part of the DRAIN mode.

RED LED LIGHT (630nm)


Photobiomodulation that supports cellular activity and collagen synthesis, helping to smooth the appearance of fine lines and wrinkles and improve firmness. It is part of the FIRM mode.

GREEN LED LIGHT (530NM)

Photobiomodulation that acts on melanin production, helping to balance skin tone and reduce dark spots, while calming the appearance of redness. It is part of the DRAIN and GLOW modes.

NEAR INFRARED LIGHT (850NM)

Long-wave light, invisible to the human eye, reaches the deepest layers of the skin, where it supports regeneration and enhances the effect of other technologies. Integrated into GLOW mode.



THREE MODES. PROGRESSIVE TREATMENT.

DRAIN
DEPUFFS. DEFINES. PREPS.
EMS + GREEN LED + VIBRATION.

- Helps reduce the appearance of swelling and redness.
- It helps to decongest and define the jawline and facial contour.
- It leaves the skin looking fresher and more revitalized.

FIRM
FIRMS. SCULPTS. SMOOTHS.
BIPOLAR RF + RED LED + VIBRATION.

- It helps improve the appearance of skin firmness and elasticity.
- It supports the skin's natural collagen.
- It helps to smooth the appearance of fine lines and wrinkles.

GLOW
ILLUMINATES. EVENS. REVITALIZES.
GREEN LED + NIR + VIBRATION.

- Enhances skin radiance.
- It helps to even out the appearance of the skin tone.
- It promotes a healthier and more radiant-looking skin.

EFFECTIVENESS

- **36% fewer visible wrinkles** in just 4 weeks[1]
- Up to **76% more newly formed collagen** at 3 months[2]
- **Visibly denser and firmer skin** in just 8 weeks[3]
- **Visibly improved tone, brightness and appearance of dark spots** in 8 weeks[4]
- Up to **+102% elastin** in 3 months[5]
- **100% of users noticed an improvement in wrinkles** after 8 weeks[3]

- [1] Randomized, double-blind, placebo-controlled, split-face, in vivo clinical study; 76 participants, 2 sessions/week for 4 weeks with red (633 nm) and near-infrared (830 nm) light. Wrinkles and elasticity were measured instrumentally (profilometry and cutometer). Individual results may vary.
- [2] Independent in vivo clinical and histological study; 6 volunteers, 6 radiofrequency sessions over 3 months. Collagen quantified in skin biopsies using histometric and immunohistochemical techniques. Individual results may vary.
- [3] Open-label, controlled, in-individual in vivo clinical study; 22 female volunteers (25–60 years old), home radiofrequency treatment for 8 weeks. Skin density assessed by ultrasound and wrinkles by clinical classification. Individual results may vary.
- [4] Randomized, split-face, in vivo clinical study; facial masks with 505 nm green LED versus a control without LED for 8 weeks. Melanin, luminosity, color index, and number of spots were measured instrumentally. Individual results may vary.
- [5] In vivo histological study; skin biopsies analyzed after four weekly sessions combining radiofrequency and facial muscle stimulation, evaluated at 3 months of follow-up. Collagen and elastin quantified by histology. Individual results may vary.

SAFETY AND COMFORT

- **100% safe technology:** safe and effective skin care device, CE certified.
- **Non-invasive technology:** RF, EMS, and LED/NIR light are non-invasive methods suitable for all skin types. Gentle heat and stimulation are kept at comfortable levels, and the light does not emit UVA rays.
- **High-quality materials:** skin-friendly contact surfaces that glide smoothly over face and neck without pulling or marking.
- **No recovery time:** unlike laser or peels, you can resume your daily activities right after each session.

UNICQUANTUM-V GuaSha™ 5.0 USAGE SCHEME

SHORT SIDE OF SCULPTING

It is used in areas where a precise line needs to be drawn following the bone because its shape fits better to those curves to mark and define the facial contour.



ON/OFF AND MODE SELECTION BUTTON

T**op button:** Press and hold for 1 second to turn on the device.
Bottom button: changes mode, in order DRAIN → FIRM → GLOW.



LONG SIDE MASSAGE

Its shape distributes pressure evenly during long, continuous glides. It's the general "drag" side, designed for drainage and relaxation rather than definition.

LED INDICATOR

LED display showing the usage of each mode (LED, DRAIN·EMS, FIRM·RF, GLOW·NIR) and power level bars for each.

USB CHARGING PORT

100% wireless. It recharges via USB-C and takes approximately 2.5 hours for a full charge. A single charge allows for multiple sessions.

HOW TO USE




STEP 1

Before starting the treatment, thoroughly cleanse and dry your face and neck to remove any traces of makeup, oil, or impurities. Remove any jewelry (earrings, necklaces, piercings) before using the device, as it could interfere with the radiofrequency and EMS contact.



STEP 2

Apply a generous, even layer of UNIC**LIFTING+** conductive serum to the entire treatment area (face and neck). The serum allows the device to glide smoothly and ensures that the radiofrequency, EMS, and LED/NIR light are properly delivered to the skin.



STEP 3

Turn on the device by pressing and holding the top button for 1 second. Press the bottom button to select the treatment mode: DRAIN, FIRM, or GLOW. Each mode combines different technologies and corresponds to a phase of the facial ritual.



STEP 4

Place the device in full contact with the skin and glide it following the movements indicated for the selected mode. Keep the device moving throughout the session and repeat each stroke 10 to 15 times per side.

After completing the treatment, massage the remaining serum into your skin until fully absorbed and continue with your usual beauty routine. Do not use the device on dry skin, as this may cause friction and reduce the treatment's effectiveness.

Complete session: 6 minutes total (2–3 minutes per mode), always following the order **DRAIN → FIRM → GLOW** for a progressive treatment: first, the skin is decongested and prepared, then firmness and contouring are addressed, and finally, radiance is enhanced. Recommended frequency: daily use, especially in the morning, as part of your facial care routine.

Mode 1 · DRAIN (2–3 min) · Decongests, refreshes and revitalizes



Neck
Begin on both sides of the neck, sliding from under the ears towards the collarbones, using the long side of the massage.



Nose
Glide outwards and slightly downwards, from the nasolabial fold towards the cheek and ear area.



Cheeks
Gently slide from the inner part of the cheek, next to the nose, outwards and slightly downwards towards the ears.



Jaw
Start in the center of the chin and glide along the jawline towards the ears with the sculpting side; then go down the neck towards the collarbones.

Gentle pressure | Towards the drainage channels | Repeat 10-15 times per side

Mode 2 · FIRM (2–3 min) · Reaffirms, sculpts, and defines



Neck
Start at the base of the neck and slide upwards towards the jaw, using the long side for massaging.



Jaw
Start at the center of the chin and slide upwards and outwards along the jawline towards the ears.



Cheeks
Start next to the nose and slide diagonally, upwards and outwards, across the cheeks towards the temples.



Eyes and forehead
Start at the inner corner of the eye and, with very gentle pressure, slide outwards and upwards towards the temples; then from the eyebrows towards the hairline.

Light-medium pressure | Upward and outward movements | Repeat 10-15 times per side

Mode 3 · GLOW (2–3 min) · Illuminates, softens and revitalizes



Neck
Start at the base of the neck and slide upwards towards the jaw with smooth, continuous movements.



Jaw
Start in the center of the chin and glide outwards along the jawline towards the ears with the sculpting side.



Forehead
Start between the eyebrows and glide upwards and outwards towards the hairline and temples.



Cheekbones
Start near the center of your cheek and gently glide outwards towards your cheekbone and temple with the sculpting side.

Gentle pressure | Smooth and fluid movements | Repeat 10-15 times per side

UNICQUANTUM-V GuaSha™ 5.0 ACCESSORIES

UNICQUANTUM-V GuaSha™ 5.0 comes in a compact and elegant case that includes everything needed to start treatment comfortably, safely and efficiently.


Device


Travel case


USB-C charging cable


Express user manual

LOAD

Before first use, fully charge the device. Completely wireless, it recharges via the included USB-C cable in approximately 2.5 hours, and a full charge allows for multiple sessions. During charging, the indicator light flashes; once charging is complete, it remains on steadily. When the battery is low, the device emits a beep and automatically shuts off. Ensure the connector is dry before charging and avoid contact with water.

CARE AND CONSERVATION

After each use, clean the device with a soft cloth dampened with water and, if necessary, a small amount of mild soap; then dry it with a clean cloth. Avoid using alcohol, acetone, or solvents on the device, as these could damage its finish. Once clean and dry, store it in its velvet case to protect it from dust and bumps and keep it in perfect condition.

TECHNICAL SPECIFICATIONS

- **Power**: 5.5 W
- **Intensity**: 4 adjustable levels
- **Wavelengths**: Red LED 630 nm · Green LED 530 nm · NIR 850 nm
- **Number of LEDs**: 8 RGB + 6 IR
- **Radiofrequency**: 500 kHz · up to 4.5 W
- **RF temperature**: controlled in 4 levels (up to 40–42 °C)
- **EMS**: frequencies of 10 / 20 / 52 Hz · 4 intensity levels
- **Battery**: Rechargeable lithium 3.7V · 1000mAh
- **Charging**: USB (5V) · ~2.5 hours for a full charge
- **Not water resistant**



PRECAUTIONS

Do not use in case of:

- Pacemaker, defibrillator, or other implanted electronic device.
- Pregnancy or breastfeeding (safety of use during these stages has not been evaluated).
- Epilepsy or a history of seizures.
- Coagulation disorders or anticoagulant treatment.
- Active cancer, tumors, or undiagnosed skin lesions.
- Wounds, infection, active herpes, eczema, or recent scarring in the area.
- Recent cosmetic procedure (laser, peeling, surgery or injectables): wait for the skin to recover and respect the times indicated by the doctor.

Consult a medical professional if you experience:

- Heart disease, thyroid disease, or autoimmune disease.
- Rosacea, couperose, melasma or tendency to hyperpigmentation.
- Hypersensitivity or reduced sensitivity to heat.
- Fillers, Botox/neuromodulators or lifting threads (see FAQ).
- Large metal implants, plates or screws in the treatment area.

Avoid applying it to:

- Eye contour and eyelids.
- Central and anterior area of the neck (thyroid) and cervical arteries.
- Chest and décolletage.
- Raised moles, warts, or prominent capillaries.

Usage recommendations:

- Always apply with a water-based conductive serum (UNICLIFTING+), never on dry skin.
- Keep the device moving, without stopping it on the same point.
- Start at the lowest intensity. Gentle heat from radiofrequency, slight tingling from EMS, and vibration are normal; reduce the intensity or stop if you experience any discomfort.
- Apply sunscreen in the following days.
- Frequency: Suitable for daily use. The FIRM (radiofrequency) mode can be used daily at low intensity; at higher intensities, it is recommended up to 3 times per week.

Equipment and conservation:

- Not water resistant: do not submerge, do not rinse or wet the charging port.
- Charge only with the included USB-C cable, with the connector dry.
- Do not disassemble the device.
- Store in a cool, dry place (max. 42°C), out of reach of children.

FREQUENTLY ASKED QUESTIONS

What is UNICQUANTUM-V™ GuaSha 5.0 and how does it differ from a traditional Gua Sha?
It's the technological evolution of Gua Sha: it takes the traditional massage and enhances it with five technologies (Bipolar Radiofrequency, EMS, and Red, Green, and Infrared LED light). Where traditional Gua Sha only works on the surface, this also works deep down, improving firmness, contour, tone, and radiance, in a facial ritual of just 6 minutes.

Why choose it over other competing facial devices?
Because it accomplishes in a single device what others require several to achieve. Unlike devices that combine fewer technologies, the UNICQUANTUM-V™ GuaSha 5.0 brings together 5 technologies, with 3 wavelengths, EMS, and adjustable radiofrequency, to simultaneously treat firmness, contour, tone, and radiance. Furthermore, its ergonomic Gua Sha shape adapts to different areas of the face and helps to precisely define the contour. A more complete treatment, in less time, and in a single device.

When will I see results?
From the first use: skin appears less congested, more radiant, and with an immediate lifting effect. With continued use, the results build, becoming firmer and smoothing the appearance of wrinkles week after week.

How often should I use it?
The complete ritual takes only 6 minutes (2-3 minutes per mode) and you can enjoy it daily. Its radiofrequency is of gentle intensity, designed for home use: at low intensity, the FIRM mode is perfect for everyday use; if you prefer higher intensities, we recommend spacing it out to about 3 times per week. You decide the frequency according to your skin's needs.

Do I need to use a conductive serum/gel?
Yes, it's part of the ritual. Apply a water-based conductive serum to cleansed skin so the device glides smoothly and the technologies are delivered optimally. Its perfect partner is UNICLIFTING+: its combination of hyaluronic acid, orange stem cells, and phytoproteoglycans hydrates, firms, and provides a lifting effect. Furthermore, EMS and radiofrequency heat help the skin better absorb its active ingredients.

In what areas can I use it?
On face and neck. Avoid the eye area and eyelids, the center of the neck, and the décolletage. Always keep the device moving.

Does it hurt? Is the heat normal?
No, it's painless. In FIRM mode, you'll feel controlled radiofrequency heat (adjustable to 4 levels, up to about 42°C), the temperature that stimulates collagen without damaging the skin. With EMS, it's normal to feel a slight tingling or mild contractions: these are the micro-impulses that work the facial muscles.

Can I use it if I have cosmetic treatments (neuromodulators, fillers or lifting threads)?
Yes, you can continue using it even after undergoing cosmetic treatments, keeping in mind a few simple guidelines. After a recent procedure, wait the recovery time indicated by your doctor. Afterward, you only need to take one precaution regarding the specific treated area:
Fillers: avoid radiofrequency (FIRM mode) on that area.
Neuromodulators: avoid EMS (DRAIN mode) over that area.
Thread lifts: avoid radiofrequency and EMS on that area.
In that area, you can continue using the device turned off, like a manual Gua Sha, and use it normally on the rest of your face. If you have any questions, consult your cosmetic doctor.

Can I use it if I have dental implants or fillings?
Yes, they are compatible. Simply keep the device moving around that area, without stopping the radio frequency directly above it.

CUSTOMER SERVICE

If you have any questions or comments you would like to share with us, please write to us at the following email address: customerservices@unicskin.com.

DISCLAIMER

This product is not intended to diagnose, cure, or prevent any disease. UNICSKIN S.L.
This product does not guarantee its ability to treat any physical, mental, or skin condition. For any health concerns, consult a qualified medical professional.

WARRANTY

Except in cases of manufacturing defects, returns are not accepted for products purchased through our website (www.unicskin.com) once they have been used and opened. If the product was purchased outside of our website, please consult the warranty policy offered by the retailer, which in no case supersedes this warranty policy.

If the packaging has not been opened, UNICSKIN will accept returns within 14 days of the purchase date. The customer is responsible for all return shipping costs and risks. The device must be returned in its original packaging, along with the original proof of purchase. Please contact Customer Service at customerservices@unicskin.com or by phone at +34 912 797 377 for information on the return address.

In the event of a manufacturing defect and during the first three years from the date of purchase, the return of any of our Techno-Beauty products may result in a replacement with the same product (in which case the return costs will be borne by UNICSKIN), a refund, or a price reduction if the defect detected is of little importance or the value of the item at the time of the claim does not justify its replacement due to the time elapsed.

In any case, it will be necessary to keep and present the receipt or invoice for the purchased product, clearly showing the purchase date. Furthermore, UNICSKIN reserves the right to carry out the necessary checks on the purchased product, requesting photographs or evidence of the alleged nonconformity or the return of the defective product. Initially, the customer will be responsible for the shipping costs, which will be refunded once UNICSKIN has verified that the detected fault corresponds to a manufacturing defect.

Products damaged due to storage conditions, maintenance, or any other type of handling by the customer are excluded from this warranty policy. The product complies with the legal and regulatory standards applicable in the European Community at the time of manufacture and distribution.

UNICSKIN

YOUR EFFECTIVE SOLUTION

GLOBAL EXPERTS IN TECH BEAUTY & HOME BEAUTY

EFFECTIVE RESULTS PROVEN BY SCIENCE

SUSTAINABILITY & FORMULAS +92% OF NATURAL ORIGIN

INNOVATIVE LUXURY WITHIN YOUR REACH
